## TECHNOCOL 6000

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### GB User instructions

#### Material

TECHNOCOL 6000 is a rapid-curing resin on a methyl methacrylate basis supplied as powder and liquid.

#### Indication

**Claw treatment:**  
Immobilisation of claws during inflammatory diseases. Immobilization of claw tips.

#### Treatment of bone fractures in large and small animals:

Extra-cutaneous osteosynthesis splinting after Becker, transversal fixation by means of Kirschner drill wires, radius and ulna fractures, tibia-fibula fractures, pedal bone fractures. strengthening of dressings.

#### Properties

TECHNOCOL 6000 is lor 3 phase processing (lor casting, lilling and modelling). When it comes to modelling TECHNOCOL6000 does not become sticky and can be modelled to any form desired without the use of aids. The curing time is about 4 - 6 min, according to mixing ratio and product temperature. Once cured TECHNOCOL 6000 has a non-sticky surface. TECHNOCOL 6000 is resistant to acids, alkalis and external impurities and can be dissolved in acetone. In a fluid and doughy state it attacks rubber. Polyethylene and Hostaphan films are not attacked. and can therefore be used to insulate and cover the dough. After completion of polymerisation TECHNOCOL 6000 becomes hard and abrasion-resistant allowing for mechanic machining (grinding, sawing, drilling, milling etc.)

#### Processing

**Preparation**  
The surface must be roughened (for example, using a hoof rasp, hoof knife or abrasive disc). The parts getting in touch with TECHNOCOL 6000 (claw, wood en block) must be free of dust and grease and completely dry (we recommend e.g. spirit or alcohol). To achieve reasonable curing times of approx.4 - 6 min., it may prove practical to cool the mixing vessel and material in summer or to heat it (to about 20–23°C) in winter. At very low outside/product temperatures

it is recommended adding TECHNOCOL ACTIVATOR to the mixture.

#### Processing with ACTIVATOR

TECHNOCOL ACTIVATOR is a co- initiator mixture developed to shorten polymerisation of TECHNOCOL 6000 at low outside/product temperatures. When outside temperatures are over + t 0°C mix TECHNOCOL 6000/8000 liquid and TECHNOCOL ACTIVATOR in a ratio of 14:1. At temperatures below +10°C the mixing ratio is 7:1. This corresponds to 500 ml of TECHNOCOL 6000 and 35 ml respectively 70 ml of TECHNOCOL ACTIVATOR.

#### Mixing TECHNOCOL 6000 - Mixing ratio of powderfluid about 2:1

Shortly before use pour fluid in a mixing beaker (see graduation). Add powder until the mass reaches a still pourable consistency. The mixing ratio can be varied according to individual requirements: Add more powder to obtain a more solid consistency and more fluid for a thinner mixture. Directly after mixing TECHNOCOL 6000 is still fluid enough to be poured. Then the resin turns into a dough-like consistency and becomes ready for filling, and a bit later, for modelling.

#### Application

Spread a thick layer of TECHNOCOL 6000 on the claw sole and wood block. Having pressed the wood block hard onto the claw sole, the surplus resin has to be modelled over the sides of the claw and the wood block. The circumferential groove of the wood block has to be filled properly with resin to provide good stability. After about 4-6 minutes according to mixing ratio and product temperature TECHNOCOL 6000/8000 will cure completely and once cured will withstand high loads.

#### Curing

The curing process begins when powder and liquid are mixed. Curing time can be shortened by heating or prolonged by cooling, but not interrupted. Curing brings forth heat. The resulting polymerisation temperature depends on the quantity of TECHNOCOL 6000 added. To avoid possible tissue damages, the resin can be sprinkled with cold water until curing is completed. The synthetic resin can be cut with a hot spatula (caution - flammable!). Openings produced in this way can be easily and permanently closed with freshly-mixed resin. When the healing process has been completed, the plastic can be mechanically removed with the use of appropriate tools (pliers, hammer).

#### Application

**Claw treatment:**  
Claws are being immobilized during inflammatory diseases by fitting a wood block to the unaffected neighbouring claw using TECHNOCOL 6000 as connector.

#### Bone Treatment:

To fix the fractured ends in their original position until healing is completed, the bone ends are connected by percutaneous bone screws

(Becker) or Kirschner drill wires and individually shaped bridge splints made of TECHNOCOL 6000. Advantages: The low weight is favorable for small domestic animals and furthermore, there is the possibility of repositioning the fragments as long as the callous hasn't hardened yet (by separating and renewed bonding over the length necessary to ensure bone healing).

#### Danger/Safety advices

Please follow the instructions on the product packaging and safety data sheets.

#### Storage advice

Do not store the material above 25°C. Store in a tightly sealed container. Do not use the material after the expiry date. Keep out of the reach of children.

#### Field of responsibility

Our usage instructions in spoken and written form, resulting from test carried out by us, are made to the best of our knowledge. However, they remain without obligation, also with reference to any trade mark rights of a third party. This does not free you from checking the products delivered by yourself that they are fit for the intended application and use. We accept no liability with regard to the application, use and processing of the product, since such use is beyond our control remaining exclusively in your field of responsibility. Of course we guarantee for the excellent quality of our products in conformity to our General Terms and Conditions of Delivery.

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